Influenza A+B

(Cassette)

REF AXP-242001-10K-W
AXP-242001-10K-O

IVD 

| | | |
|-----------|----------------------------|----|
| DE | Gebrauchsanweisung | 2 |
| EN | Instructions for use | 5 |
| FR | Instructions d'utilisation | 8 |
| IT | Istruzioni per l'uso | 11 |
| SI | Navodila za uporabo | 14 |
| CZ | Návod k použití | 17 |

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14. Referenzen

- Murphy, B.R., and R.G. Webster, 1996, Orthomyxoviruses, pp.1397-1445. In: Fields, Virology, 3rd edition, B.N. Fields, D.M. Knipe, P.M. Howley, et al. (eds.), Lippincott-Raven, Philadelphia.
- Biosafety in Microbiological and Biomedical Laboratories, 4th Edition. U.S. Department of Health and Human Services, CDC, NIH, Washington, DC (1999).
- CDC website: <http://www.cdc.gov/flu/>

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1. Intended Use

The AXpress® Influenza A+B test is a rapid chromatographic immunoassay for the qualitative detection of Influenza type A and B antigens (nucleoproteins) extracted from respiratory specimens. The test is designed to detect nucleoproteins of core subtypes, independent of H and N surface proteins of Influenza A (including H1N1 and H3N2). For professional use only.

2. Introduction and Clinical Significances

Influenza is a highly contagious viral infection of the upper respiratory tract which is characterised by antigen variability, seasonality and the impact on the general population. Of the two main types (A and B) of influenza viruses, Influenza A subtypes are differentiated by antigen variability of the surface glycoproteins (hemagglutinin and neuraminidase). Influenza A virus is the most prevalent and is associated with the most serious epidemics. Influenza can cause severe complications such as bronchitis or pneumonia, particularly in children, elderly people or people with chronic respiratory disease. Yet, it most commonly occurs as a mild viral infection transmitted by respiratory secretions through sneezing or coughing. There are many other viral infections that can mimic influenza symptoms, making laboratory tests necessary to distinguish it from other acute respiratory infections. New effective antivirals have been available since the late 1990s, but these drugs are only effective if administered early (before 48 hrs after the onset of the illness). Virus isolation is still considered to be the gold standard method for influenza diagnosis, with a sensitivity of nearly 100% after 3 days. However, patient health care and economic costs could be greatly improved by the use of rapid, specific and sensitive tests for antigen detection in order to enable a successful antiviral treatment.

3. Principle of the Test

The AXpress® Influenza A+B test is a qualitative, lateral flow immunoassay for the detection of Influenza A or Influenza B antigens (nucleo-proteins) in nasal swab specimens. In this test, antibodies specific to antigens of Influenza A and Influenza B strains are precoated in the test line region of the test cassette. During testing, the extracted antigen solution reacts with antibodies to the respective Influenza type conjugated with particles. The mixture migrates along the membrane, reacts with antibodies to the corresponding Influenza type and generates a coloured line in the test line region. The presence of this coloured line in the test line region indicates a positive test result, while its absence indicates a negative test result. To serve as a procedural control, a coloured line should always appear in the control line region, indicating that the proper volume of specimen has been added and membrane wicking has occurred.

4. Reagents and Materials Supplied

- 10 AXpress® Influenza A+B Test Cassettes incl. single use droppers
- 1 extraction buffer bottle
- 10 extraction tubes
- Provided additional material according to 39/42/EEC: 10 sterilised swabs for sample collection CE 0197

Dalian Rongbang Medical Healthy Devices Co. Ltd
Maoyingzi Hamlet Dalianwan Town
Gangjingzi District 116113, P.R.C.
(authorized EU representative: Yaotong S.L.
Ausias Marc, 92-98, esc. C, entre. 3ª
Barcelona, 08013, Spain).

- 1 Package insert

5. Additional Required Materials

- Timer

6. Storage & Stability

The AXpress® Influenza A+B test should be stored as packaged in the sealed foil pouch either at room temperature or refrigerated (2-30°C). The test is stable until the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. Do not freeze. Do not use beyond the expiration date.

7. Warnings and Precautions

- For professional in-vitro diagnostic use only.
- For single use only.
- Do not eat, drink or smoke in the area where the specimens or test kits are handled.
- All work steps related to the use of the test should be performed in compliance with the requirements of Good Laboratory Practice (GLP).
- Do not use after the expiration date, i.e. no guarantee of quality can be given after the expiry date or in the case of improper storage.
- Bring specimens and test components to room temperature (15-30°C) before testing.
- Do not use the test if the pouch is damaged.
- In order to avoid contamination, do not touch the reaction (result) window, i. e. do not touch the nitrocellulose membrane with your fingers.
- Do not add specimen samples to the result window.
- Please use a new specimen collection tube for each specimen in order to avoid cross-contamination.
- Do not mix reagents from different lots, i.e. reagents from different test kits should not be used.
- Do not use extraction buffer contaminated with bacteria or mould.
- Handle all specimens as if they contain infectious agents. Observe established precautions against microbiological hazards throughout testing and follow standard procedures for proper disposal of specimens. Wear protective clothing such as laboratory coats, disposable gloves and eye protection when specimens are being tested. The used test should be discarded according to local regulations.
- Humidity and high temperature can adversely affect test results.
- Perform the test immediately after removing it from the foil pouch.
- Extraction buffer contains an antimicrobial agent that does not present any danger as long as the normal laboratory safety precautions are observed.
- Avoid splashing or aerosol formation.
- Clean up spills thoroughly using an appropriate disinfectant.

8. Specimen Collection and Preparation

Nasal Swab

- To collect a nasal swab specimen, insert the sterile swab into the nostril producing the most secretion under visual inspection.
- Using gentle rotation, push the swab until resistance is met at the level of the turbinate (less than one inch into the nostril).
- Rotate the swab a few times against the nasal wall.




Sample Storage & Transport

The sample should be tested as soon as possible after collection in order to obtain the most accurate test results.

9. Procedure of the Test

Bring the AXpress® Influenza A+B test, specimens and controls to room temperature (15-30°C) prior to testing.

- Add 15 drops of extraction buffer to the extraction tube.
- Immediately insert the swab with collected nasal swab specimen into the tube, swirl the swab vigorously swirl the swab vigorously for at least 1 minute in order to mix the reagents. Squeeze as much liquid as possible from the swab by rotating and pressing its fiber tip against the wall of the tube. Discard the swab. Mix the contents of the tube by gentle swirling. The mixture is ready to be tested.
- Remove the AXpress® Influenza A+B test cassette from the foil pouch and place it on a flat and dry surface.
- Holding the sample dropper above the test cassette, add 5 drops of the mixed specimen to each specimen well. Wait until each drop is absorbed before adding the next drop. If there is no colour migration across the result window within approximately 30 seconds, add additional drop(s). As the test begins to work, you will see red colour move across the result window in the centre of the test cassette.
- Interpret test results after 10 minutes. Do not interpret test results after 15 minutes.

Note: The above interpretation time is based on reading test results at room temperature of 15 to 30°C. If your room temperature is significantly lower than 15°C, then the interpretation time should be properly increased.

10. Interpretation of the Results

Positive

Generally, the higher the analyte level in the sample, the stronger the colour of the line in the test line region (T). If the analyte level in the sample is close to the cut-off value (the detection limit of the test), the colour of the line in test line region (T) will be very faint.

Note:

A positive result will not change delete after it has been determined at 15 minutes. However, in order to prevent incorrect results, test results should not be interpreted after 15 minutes.

Positive for Influenza A

One red line develops in each control line region (C) and one red line develops in the test line region (T) for Influenza A.

Positive for Influenza B

One red line develops in each control line region (C) and one red line develops in the test line region (T) for Influenza B.

Positive for Influenza A + B

In addition to the control lines (C), a red line develops in each test line region (T) for Influenza A and Influenza B.


Negative

One red line develops in each control line region (C). No line develops in either test line region (T).

Invalid

The absence of the control lines (C) indicates a failure in the test procedure. The sample should be retested.

Note: Insufficient specimen volume, incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your distributor.

11. Quality Control

Control standards are not supplied with this test kit; however, it is recommended that external controls be tested as Good Laboratory testing practice.

12. Limitations

Although the AXpress® Influenza A+B test is accurate in detecting Influenza A/B virus, false results can occur. Other clinically available tests are required if questionable results are obtained.

The AXpress® Influenza A+B test is a qualitative assay. The amount of Influenza A or B present in the specimen cannot be estimated by the assay. The assay results distinguish positive from negative samples only.

As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of a single test, but should only be made by the physician after all clinical and laboratory findings have been evaluated.

13. Performance Characteristics

Sensitivity and Specificity

The AXpress® Influenza A+B test was compared against RT-PCR for its relative sensitivity and specificity: The RT-PCR relative sensitivity and specificity of the AXpress® Influenza A+B test is 84.4%/94.7% for A strains and 84.3%/95.5% for B strains.

Sensitivity (Analytical Sensitivity)

Virus titration in embryonated chicken eggs was used to determine analytical sensitivity.

| Virus strain | Type | Subtype | CEID/Test |
|-----------------|-----------|---------|------------------------|
| Texas 1/77 | A | H3N2 | 2.9 x 10 ³ |
| HongKong/8/68 | A | H3N2 | 6.1 x 10 ¹ |
| Shanghai/11 | A | H3N2 | 6.1 x 10 ¹ |
| Port Chalmers | A | H3N2 | 1.25 x 10 ² |
| New Jersey/8/76 | A (Swine) | H1N1 | 2.5 x 10 ² |
| NWS/33 | A | H1N1 | 1.24 x 10 ³ |
| Mal/302/54 | A | H1N1 | 6.73 x 10 ³ |
| HongKong/5/72 | B | N/A | 3.6 x 10 ¹ |
| Taiwan | B | N/A | 3.6 x 10 ² |
| California | B | N/A | 3.1 x 10 ³ |
| Stockholm | B | N/A | 3.6 x 10 ⁵ |

In addition, the AXpress® Influenza A+B tests were shown to detect all the examined strains.

H2N2 – A/Hokkaido 17/01 (duck)
 H3N8 – A/Mongolia/3/03 (duck)
 H5N1 – A/HongKong/156/97
 H5N1 – A/Yamaguchi/7/04 (chicken)
 H5N1 – A/Thailand/MK2/04
 H6H2 – A/Massachusetts/3740/65 (turkey)
 H7N7 – A/Massachusetts/1/80 (seal)
 H9N2 – A/Wisconsin/66 (turkey)

14. References

- ① Murphy, B.R., and R.G. Webster, 1996, Orthomyxoviruses, pp.1397-1445. In: Fields, Virology, 3rd edition, B.N. Fields, D.M. Knipe, P.M. Howley, et al. (eds.), Lippincott-Raven, Philadelphia.
- ② Biosafety in Microbiological and Biomedical Laboratories, 4th Edition. U.S. Department of Health and Human Services, CDC, NIH, Washington, DC (1999).
- ③ CDC website: <http://www.cdc.gov/flu/>

Abbreviations/Symbols

| Symbol | Meaning |
|---|---------------------------------------|
|  | Catalogue number |
|  | Contents of kit |
|  | For in vitro diagnostic use |
|  | Manufacturer |
|  | Storage temperature |
|  | Use before |
|  | Lot number |
|  | Consult the package insert before use |
|  | Single use |
|  | CE marking |
|  | Sufficient for <n> tests |

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nal von minden GmbH
Axon Lab AG